

## When to use Parafricta® bedding

For fragile skin conditions such as

- at high risk of pressure ulcers
- early stage pressure ulcers
- redness and friction lesions due to agitation
- treated burns or skin grafts
- epidermolysis bullosa
- allergic or autoimmune skin reactions
- end of life skin breakdown

To aid in turning over in bed

- in neurodegenerative conditions such as Parkinson's
- in arthritic conditions

## Full Instructions for use and product dimensions

Please see User Guide



## Washing Instructions

Please see Care Label or User Guide

[www.parafricta.com/user-guides-27-w.asp](http://www.parafricta.com/user-guides-27-w.asp)


**Parafricta®**

APA Parafricta Ltd

01869 226566

[info@parafricta.com](mailto:info@parafricta.com) [sales@parafricta.com](mailto:sales@parafricta.com)

[www.parafricta.com](http://www.parafricta.com)



**Parafricta®**

**BEDDING**

***Non-Stop Relief from Friction  
and Easier Turning in Bed***

Parafricta® Bedding employs a unique low friction fabric over the main area and is designed to

- **protect fragile skin** in contact with the support surface against the damage caused by friction and associated shear
- improve **comfort** of users
- and/or aid users with limited mobility to **self-reposition** (turn) whilst in bed.

**There is a range of Parafricta® sheets and pillowcases for use in different settings.**

All bedding can be washed at up to 95°C (same as regular bedding), tumble dried and gently ironed. Gives cool, silk-like comfort against the skin.



\* The picture on the cover shows a blue border for illustrative purposes only. This area is made of polycotton to tuck in. The products are all white.

## Which product should I use?

Consider these factors

- (1) the reason for using the sheet (a) skin protection (b) aid to self-repositioning
- (2) the type of bed frame i.e. flat or tilting/profiling
- (3) the size of the mattress (a) single (b) double/king
- (4) the depth of the mattress

## Select the right product

	Product Code	Used for	Type of bed or pillow	Mattress Size
Single fitted sheet	PFFSS-001	Skin protection	Only flat bed	Standard single
Single flat sheet	PFFSS-003	Skin protection	Flat or tilted/profiling bed	Single, including deep*
Single, semi-fitted sheet (Mobility Assist)	PFFMAS-001	Self-repositioning	Flat or tilted/profiling bed	Standard single
Double flat sheet (Dual Purpose)	PFDDS-001	Self-repositioning or Skin protection	Flat or tilted/profiling bed	Double or King, including deep
Pillowcase	PFFPC-001	Skin protection	UK rectangular pillow	n/a

\*Flat sheets may be used over dynamic mattresses but should not be tucked in tightly